

Product Overview

The percutaneous insertion kit is designed to facilitate percutaneous insertion of a cannula. The kit includes the following components:


- Five dilators: 8 Fr / 12 Fr / 16 Fr / **20 Fr** / 24 Fr
- Number 11 scalpel
- 18 ga. insertion needle
- 5 mL syringe
- 0.038 inch guidewire

5 units per case

Arterial Insertion Kit

PIKA5


100 cm guidewire



Venous Insertion Kit

PIKV5

210 cm guidewire



For more information:

Contact your sales rep | Call 1-888-943-2783 | Email customerservice@enablecv.com

CAUTION: Federal (United States) law restricts this device to sale by or on the order of a physician. See instructions for use for full prescribing information, including indications, contraindications, warnings, precautions, and adverse events.

enableCV and the stylized enableCV logo are trademarks of enableCV, Inc. and its affiliates. All other trademarks are the property of their respective owners.

© 2025 enableCV, Inc. All rights reserved. ECV-SM-0024 01