


Product Overview

Innovative design for ease of use and uncompromised flow¹

Indication

The QuickDraw Plus venous cannula serves to drain non-oxygenated blood from the venae cavae or right atrium during cardiopulmonary bypass.

Features of the QuickDraw Plus Femoral Venous Cannula



1. enableCV bench data on file, may not be indicative of clinical performance.

Ordering Information

QuickDraw Plus Femoral Venous Cannula

Model # **QDP25**

Package contents

- 25 Fr (8.3 mm) cannula (65 cm length)
- Introducer with alignment marker
- 3/8 in (9.5 mm) barbed connector hub


Model # **QDP22**

Package contents

- 22 Fr (7.3 mm) cannula (65 cm length)
- Introducer with alignment marker
- 3/8 in (9.5 mm) barbed connector hub

Note: Insertion components are not included.

PIKV5 Percutaneous Insertion Kit is recommended.




**Mean value derived from In-Vitro testing performed with water at ambient temperature. Due to increased viscosity of blood, actual pressure gradients encountered in a clinical situation may vary from these shown.*

Percutaneous Insertion Kit Venous Insertion

Model # **PIKV5**

Package contents

- Five dilators: 8 Fr / 12 Fr / 16 Fr / 20 Fr / 24 Fr
- Number 11 scalpel
- 18 ga. insertion needle
- 5 mL syringe
- 0.038 inch guidewire, length 210 cm



For more information:

Contact your sales rep | Call 1-888-943-2783 | Email customerservice@enablecv.com

CAUTION: Federal (United States) law restricts this device to sale by or on the order of a physician.



See instructions for use electronically at eIFU.enableCV.com for full prescribing information, including indications, contraindications, warnings, precautions, and adverse events.

enableCV, the enableCV stylized logo, and QuickDraw are trademarks of enableCV, LLC and its affiliates. All other trademarks are the property of their respective owners.

© 2026 enableCV, LLC. All rights reserved. ECV-SM-0026 04



enableCV • enablecv.com
6967 S. River Gate Drive, Suite 104, Midvale, UT 84047